

Protocol EMULSIFLEX C5

▶ Suggested lysis buffer (**mandatory**):

50 mM Tris-HCl pH8.0
50-100 mM NaCl
1 mM EDTA
10 mM MgCl₂
DNase @ 0.1mg/mL
Lysozyme @ 0.1mg/mL
Protease Inhibitors

- ▶ Incubate the bacterial cell suspension for 30min @ 4°C with gentle shaking.
- ▶ Wash the Emulsiflex “sample cylinder body” (SCB) with 200 mL dH₂O, to remove Ethanol 20%.
- ▶ Switch on the pressure knobs (**PK**) up to 10-15 psi. Output pressure should be between 10000-15000 psi for *E. coli* cells (same OD as for the French Press).
- ▶ Equilibrate the SCB with your lysis buffer (50 mL). Pressure will barely increase with water/buffer.
- ▶ Transfer your sample into the SCB (7 to 280 mL samples).
- ▶ 2-3 cycles of lysis at 10000-15000 psi. Let the lysate going through the tubing by adding 10-20 mL of lysis buffer.
- ▶ Wash with 250 mL dH₂O (pressure ON) and clean up the inside wall with water.
- ▶ Wash with 50 mL NaOH 0.5 M (pressure OFF).
- ▶ Wash with 250 mL dH₂O (pressure OFF).
- ▶ Wash with 100 mL EtOH 20% (pressure OFF). Leave about 25 mL in the system tubing.
- ▶ Empty the waste...

REMs:

- Please fill up the various cleaning solutions when they are near empty. Many thanks.
- We have an adaptor to mount syringes (20 and 50 mL) on the Emulsiflex. Useful for working with small volume samples.
- Please put DNase as DNA clumps will probably block the machine for a long time...

Any question: edurand@imm.cnrs.fr